

# Univerzalni merač temperature

Medicinski infracrveni termometar

Merenje u  
1 sekundi

Merenje  
telesne temperature

Merenje  
temperature objekta

Merenje  
temperature mleka

Merenje  
temperature vode



А К А Д Е М И Ј А  
ФИЛИПОВИЋ

ul. Kneza Lazara L2 +381 63 824 99 99  
lokal 15,16 Jagodina +381 35 82 00 999  
akademijafilipovic@gmail.com



normalna

niska

visoka



≤ 37.5°C



37.6°C - 38°C



≥ 38.1°C

## Merenje u jednoj sekundi



# Bez kontaktni infracrveni termometar

Telesna temperatura: 32.0°C - 42.9°C

Temperatura objekta: 0°C - 50°C

Detekcija u sekundi

Upozorenje na groznicu

Greška 0.2°C



Malih dimenzija,  
lak za držanje

Oblikovano dugme,  
pogodno i komforno  
za korišćenje  
100000 puta

